

# WED-180V

Full Digital Ultrasonic Diagnostic System

[www.weld.net](http://www.weld.net)



**WELLD**<sup>®</sup>

# WED-180V

## Veterinary Ultrasound Scanner



WED-180V is suitable to do an ultrasound examination and GA measurement for cat, dog, swine, sheep and equine.



Scanning Mode	Convex /Linear/Micro-convex
Probe Connector	2
Cine-loop	≥500 frames
Standard Configuration	6.5MHz Vet Rectal Linear Probe
Optional Configuration	3.5MHz Convex Probe
	7.5MHz Liner Probe
	5.0MHz Micro-convex Probe
	Application software
	Video Printer
	Trolley cart
	Biopsy Bracket
Image Storage	≥64 frames
Scanning Depth	≥250mm
Scanning Angle	Visible and adjustable
Display Mode	B, B+B, B+M, B+2M, M, 4B
TGCI	Total gain, 8-Segment TGC
Image Process:	Pseudo-color, image smoothen,THI, Gama correction, Histogram
Real-time Depth Measurement:	16-level adjustable
	Distance, circumference, area, volume obsterics
	(GA for equine, bovine, sheep, swine,cat, dog)
Report:	Measurement reports automatically generate
Body Marks	≥16 types
Notation	Date, time, name, patient ID, sex, age, doctor, hospital,
	full screen annotation edit
Port	PAL-D Video, RS-232, USB2.0

### Optional Probe:



3.5MHz Convex Probe



7.5MHz Linear Probe



6.5MHz Vet Rectal Probe



5.0MHz Micro-convex Probe



Distributor:

Shenzhen Well.d Medical Electronics Co., Ltd.

<http://www.welld.net> <http://www.welld.com.cn> Email:[export@welld.com.cn](mailto:export@welld.com.cn)  
 ADD: Well.D Park Qinglan 3 Rd., National Biopharmaceutical Industrial Base,  
 Pingshan New Area,Shenzhen 518118,China  
 Tel:+86-755-36900019/26073350 Fax:+86-755-36900018/26073919

