

WED-2000AV

Veterinary Ultrasound Scanner

www.wellid.net

Palm-size WED-2000AV is enjoying a great success in the global market by its innovative in design. With advantages like light weight, easy carrying, accurate imaging, simple interface and etc. Its functions like selectable two different kinds of probe, Multi-frequency offers a wide range of applications. It is really a marvelous choice for VET.



WELLD[®]

WED-2000AV

Veterinary Ultrasound Scanner



Feature

Weight 0.8kg, compact and easy carrying
 5 inch TFT LCD shows clear & stable images and reduces visual fatigue
 Built-in big-capacity lithium battery makes VET diagnosis available in no AC situation
 Advanced power management system supports long working time of battery
 Sector probe, no need distention of bladder, increase efficiency of gynecological examination
 Ergonomically design, friendly operation

Power supply

AC- adaptor or built-in rechargeable battery
 Battery can be used for more than 3 hours ' continuous working and can be charged by either AC-adaptor or car charger.



For VET

WED-2000AV is a good partner for veterinary to make regular diagnosis and calculate GA of cat, dog, swine, sheep, equine, bovine, etc.

3.5MHz Multi-Frequency Sector Probe with freeze key or water-proof function



power adapter

AC adapter with long cable makes the scanner can be used conveniently even when the battery is charged at the same time.

Technical specification

Scanning mode :	Mechanical sector scan
Probe working frequency :	2.5MHz, 3.5MHz, 5.0MHz
Image storage :	≥128 frames, no data lose when power off.
Display depth :	190 mm
Display mode :	B, B+B, B+M, M
Operation interface :	Chinese/English switchable
Image control:	Brightness, contrast ratio
Image process :	Image Smoothen, Histogram, Pseudo color
Real-time depth :	8-level adjustable
Measurement:	distance, circumference, area, volume, HR and obstetrics (GA for equine, bovine, sheep, swine, cat, dog)
Notation :	Name, age, sex, date, time, hospital, language, full screen words edit
Port :	Mouse, Video, USB2.0
Battery capability :	≥2200 mAh

Standard configuration

Main Unit
 3.5MHz multi-frequency sector probe
 Mouse
 AC-adaptor
 Lithium rechargeable battery
 Metallic suitcase

Optional

Lithium rechargeable battery
 Car charger
 power adapter
 Video printer



AISON
 INTERNATIONAL

Authorized Distributor
 AISON INTERNATIONAL
 Stigs Center 2B,
 42246, Hisings Backa, Sweden
 Tel:0046-73-5543108
 Email: info@autoclave.se
 www.autoclave.se

Shenzhen Well.d Medical Electronics Co., Ltd.

http://www.welld.net http://www.welld.com.cn Email:export@welld.com.cn
 ADD: Well.D Park Qinglan 3 Rd., National Biopharmaceutical Industrial Base,
 Pingshan New Area, Shenzhen 518118, China
 Tel:+86-755-36900019/26073350 Fax:+86-755-36900018/26073919

