

Company Milestone

2002: Company founded in Shenzhen, China

2008: Received "European Entrepreneurial Company 2008" award from FROST & SULLIVAN

2009: Received "Product Quality Leadership Award 2019" from FROST & SULLIVAN

2011: Received the Reddot 2011 Product Design Award for S20 in Essen, Germany

2013: Received "Ultrasound Market Growth Leadership Award, 2013" from FROST & SULLIVAN

2014: Received the iF product design award 2014 for S9 in Munich, Germany

2016: Received "Product Innovation in Electronic Endoscopy Market" from FROST & SULLIVAN

2017: Received the iF product design award 2017 for X5 in Munich, Germany

2017: Listed on Shenzhen Stock Exchange on April




ISO 13485

CE 0197

SonoScape

SonoScape Medical Corp.
2F, 12th Building, Shenzhen Software Park Phase II, Keji
Middle 2nd Road, Shenzhen, China
Tel: +86-755-26722890 Fax: +86-755-26722850
Email: market@sonoscape.net www.sonoscape.com

SonoScape



Comprehensive Probe Catalogue

Caring for Life through Innovation

U-TC20191024

A Wide Range of Advanced Probes

Linear Probes



12L-A

Application: Vascular, Small Parts, Breast
Musculoskeletal
Field of View: 51mm



12L-B

Application: Vascular, Small Parts, Musculoskeletal
Field of View: 38mm



9L-A

Application: Vascular, Small Parts, Musculoskeletal
Field of View: 35mm



L742

Application: Vascular, Small Parts, Musculoskeletal
Field of View: 38mm



L752

Application: Vascular, Small Parts, Musculoskeletal,
Breast
Field of View: 51mm



L741

Application: Vascular, Small Parts, Musculoskeletal
Field of View: 46mm



L743

Application: Vascular, Small Parts, Musculoskeletal
Field of View: 47mm



L746

Application: Vascular, Small Parts, Musculoskeletal
Field of View: 50mm

Micro Convex Probes



C322

Application: Abdominal Biopsy, Urology
Field of View: 68°



C611

Application: Abdomen, Cardiac, Urology, Pediatrics
Field of View: 135°



C613

Application: Abdomen, Cardiac, Urology, Pediatrics
Field of View: 90°

Convex Probes



C1-6A

Application: Abdomen, Obstetrics, Gynecology
Field of View: 70°



C1-6

Application: Abdomen, Obstetrics, Gynecology
Field of View: 64°



3C-A

Application: Abdomen, Obstetrics, Gynecology
Field of View: 60°



C344

Application: Abdomen, Obstetrics, Gynecology
Field of View: 72°



C361

Application: Abdomen, Obstetrics, Gynecology
Field of View: 72°

Volume Probes



VC2-9

Application: Abdomen, Obstetrics, Gynecology
Field of View: 63°



VC6-2


Application: Abdomen, Obstetrics, Gynecology
Field of View: 68°



VE9-5

Application: Gynecology, Obstetrics
Field of View: 182° x 120° (Volume Angle)

Phased Array Probes



S1-5

Application: Cardiac, Transcranial, Abdomen
Field of View: 90°



7P-A

Application: Cardiac (Pediatric), Neonatal Cephalic
Field of View: 90°



7P-B

Application: Cardiac (Pediatric), Neonatal Cephalic
Field of View: 90°



4P-A

Application: Cardiac, Transcranial, Abdomen
Field of View: 90°



3P-A

Application: Cardiac, Transcranial, Abdomen
Field of View: 90°



8P1

Application: Cardiac (Pediatric), Neonatal Cephalic
Field of View: 90°



2P1

Application: Cardiac, Transcranial, Abdomen
Field of View: 90°



5P1


Application: Cardiac (Pediatric), Neonatal Cephalic
Field of View: 90°



5P2

Application: Cardiac (Pediatric), Neonatal Cephalic
Field of View: 90°

Bi-plane Probes



BCL10-5

Application: Urology
Field of View: 200°/60mm



BCC9-5

Application: Urology
Field of View: 140°

Intraoperative Probes



10I2

Application: Musculoskeletal, Small Parts, Nerve, Vascular, Surgery
Field of view: 25mm



6CT-A

Application: Intraoperative, Abdomen
Field of View: 44°



12LT-A

Application: Intraoperative, Abdomen
Field of View: 33mm



6CI-A

Application: Intraoperative, Abdomen
Field of View: 44°



12LI-A

Application: Intraoperative, Abdomen
Field of View: 33mm

Endocavity Probes



6V3

Application: Gynecology, Obstetrics
Field of View: 194°



6V3A

Application: Gynecology, Obstetrics
Field of View: 147°



6V1

Application: Gynecology, Obstetrics
Field of View: 135°



6V7

Application: Gynecology, Obstetrics
Field of View: 194°



EC9-5

Application: Urology
Field of View: 147°

CW Pencil Type Probes



Application: Cardiac



Application: Cardiac, Pediatrics



Application: Neonatal cardiac , Transcranial

TEE Probes



Application: Cardiac
Field of View: 90°




Application: Pediatric Cardiac
Field of View: 90°

Laparoscopy Probe



Application: Surgery
Field of view: 39mm



*Full Range of Transducers
Provides
a Comprehensive
Clinical Solution*